


CryoHealth^o

Kloodos Wellness: Dry float

UNIQUE FEATURES

- New technology designed to simulate the effects of traditional Float-REST therapy
- No direct contact with the salt water.
- Very thin membrane that cradles the body and provides a warm sensation, similar to that experienced in a traditional flotation tank.
- Indescribable feeling of weightlessness.
- Depending on the options you choose, it also pampers you with sound waves and dry aqua massage.
- Deeply relaxing experience that is reminiscent of being in the womb.
- Reverse brain ageing

PRODUCT SPECIFICATIONS



Standard cladding

High position 1340
Low position 1110

2398

820

850

1094

Design cladding

High position 1300
Low position 1110

820

850

1090

2450

On-site connections

Wall socket 230 V, 50 Hz, 2500 W, power supply via fuse 16A.

A connection 1/2" with cold tap water for the filling is required close to the appliance.

To drain the unit: A floor drain, or an alternative way of evacuating the used water, is necessary.

For the ear-phone and option Comfort-module a microphone cable 2-wires, shielded to the head-phone jack of the on-site audio-system. Floor or wall connection.

A 2-wire cable for the connection with an on-site emergency system.

On-site safety measures

Weight of the appliance: 580 kg (including 350 litres of water) on a surface of 2 m².

With an air humidity of more than 70% a dehumidifier is necessary.

The appliance is not made for being used in aggressive air conditions (air containing salt, Sulphur or chlorine).

Minimum door width for bringing in the appliance: 70 cm. Minimum niche width: 250 cm.

Soil resistance: 340 kg/m².

9

ONLY ON DEMAND

WWW.CRYOHEALTH.BE